KATSAN PA Polyamid (PA)

Emilemeyen cerrahi iplik

TANIM

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Non-absorbable suture

DESCRIPTION

KATSAN PA sutures are synthetic non-absorbable sterile surgical sutures which are composed of the chain aliphatic polymers Polyamide 6 or 6,6. **KATSAN PA** suture is pigmented blue or black to enhance visibility. **KATSAN PA** sutures meet all requirements for synthetic non-absorbable surgical sutures specified at the European Pharmacopoeia (EP) and United States Pharmacopoeia (USP). **KATSAN PA** sutures has been presented in USP 11/0 – USP 1 gauge sizes in a variety of lengths from 10 cm to 150 cm, as non-needled or attached to stainless steel needles of varying types and sizes.

INDICATIONS

KATSAN PA sutures are indicated for general use in soft tissue approximation or ligation and ophthalmic surgery, but not for use in cardiovascular or neural tissue.

ACTIONS

KATSAN PA suture elicits a minimal acute inflammatory reaction in tissue, which is followed by gradual encapsulation of the suture by fibrous connective tissue. While nylon is not absorbed, progressive hydrolysis of the nylon in vivo may result in gradual loss of tensile strength over time. **KATSAN PA** suture provides permanent wound support.

CONTRAINDICATIONS

Due to gradual loss of tensile strength which may occur over prolonged periods in vivo, nylon suture should not be used where permanent retention of tensile strength is required.

WARNINGS

Users should be familiar with surgical procedures and techniques covering non-absorbable sutures before employing **KATSAN PA** sutures for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material is used. As with any foreign body, prolonged contact of this or any other suture with salt solutions such as those found in the urinary or biliary tracts, may result in calculus formation. Acceptable surgical practice should be followed with respect to drainage and closure of contaminated or infected wounds.

Sterility is preserved only when opened under sterile conditions.

Do not re-sterilize.

Discard opened packages and unused sutures. Do not use opened or damaged products.

PRECAUTIONS

In handling **KATSAN PA** or any other suture material, care should be taken to avoid damage from handling. Especially avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. In order to adequate knot security requires the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstance and the experience of the surgeon.

In order to prevent damaging needle points and swage areas, the needle should be grasped on its flattened area which is one-third (1/3) to one-half (1/2) of the distance from the swaged area of the needle. Reshaping needle may result in loss of strength of the needles and may result in less resistance of the needle to bending and breaking. To prevent inadvertent needle sticks, users should pay attention when handling surgical needles. Discard used needles and PA Surgical suture in "sharps" container.

Do not use expired products.

ADVERSE REACTIONS

Adverse effects associated with the use of this device include wound dehiscence, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, infection, minimal acute inflammatory tissue reaction, and transitory local irritation at the wound site.

Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood borne pathogens.

STERILITY

KATSAN PA sutures are sterilized by Ethylene Oxide.

Do not re-sterilize!

Do not use if package is opened or damaged!

Discard opened unused sutures.

STORAGE

Recommended storage conditions: Store between 5°C - 25°C,


Protect from humidity.
Do not use after expiry date!

HOW SUPPLIED

KATSAN PA sutures are available as sterile, monofilament dyed blue or black strands in sizes USP 11/0 through USP 1 (metric sizes 0.1 - 4), in a variety of lengths, with or without various needles.

Sutures are presented in one and two dozen boxes.

SYMBOLS USED ON THE LABELS

	→	CE Symbol and identification number of the notified body. The product satisfies the Essential Requirements of directive 93/42/EEC on medical devices.
	→	Manufacturer
	→	Date of manufacture (Year/Month)
	→	Use by expiration date (Year/Month)
	→	Lot Number
	→	Product Code
	→	EO : Sterilization Method - Ethylene Oxide
	→	Caution
	→	Consult instructions for use
	→	Keep away from sunlight
	→	Protect from humidity
	→	Store between : 5°C - 25°C
	→	Single use
	→	Do not re-sterilize
	→	Do not use if package is damaged

 **KATSAN KATGUD SANAYI VE TICARET A.S.**
1571 Sk. No:22 35110 Cinarlı Izmir - TÜRKİYE
Tel: +90 (232) 486 0910 Fax: +90 (232) 486 3379
Web: www.katsanas.com E-mail: info@katsanas.com