TERACEL®

Emilebilir Hemostat

Okside rejenerer selüloz

TANIM

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Oxidized regenerated cellulose

DESCRIPTION

TERACEL® Absorbable Haemostat is an absorbable and sterile woven fabric, which is generated by oxidation of regenerated cellulose. Said agent has a light yellow color and caramel-like smell. It should be stored in controlled room temperature. Its color can change in time, but this change does not have any adverse effect on **TERACEL®** Absorbable Haemostat's performance.

INDICATIONS

TERACEL® Absorbable Haemostat is used for controlling capillary bleedings and other kinds of bleedings during surgical operations. **TERACEL®** Absorbable Haemostat can be suitable for use where other kinds of haemostatic applications are not appropriate or convenient. This absorbable haemostatic agent can be used in many types of surgical methods such as; general surgery, digestive system surgery, nervous system surgery, plastic surgery, orthopaedics, gynaecology, traumatology, endoscopic operations, dental implementations and the ones related to other surgical branches.

ACTIONS

TERACEL® Absorbable Haemostat should be used singly. Absorbable haemostatic agent can be sewn, used for bandaging a bleeding wound or implicated by pressing on it. After **TERACEL®** Absorbable Haemostat is applied on bleeding surface, absorbable haemostatic agent will build up a dark gelatine mass. Thus; coagulation will occur and help bleeding to stop. When **TERACEL®** Absorbable Haemostat is used in appropriate amount required; agent will be absorbed by area of use, without any tissue reaction. Absorption of absorbable haemostatic agent by area of use depends on facts such as amount, texture and blood saturation degree of haemostatic agent used.

CONTRAINDICATIONS

TERACEL® Absorbable Haemostat is not suitable for being used as tampon or bandage. This agent can have adverse effect on callus formation and carry risk of cyst formation; so should not be used for implantation in bone damages such as breaks. **TERACEL®** Absorbable Haemostat should not be used for taking wide artery bleedings under control. Excluding blood, body fluids such as serum, do not react with **TERACEL®** Absorbable Haemostat; so it should not be used for generating haemostatic effect on surfaces, where there is no bleeding but serous leakage.

WARNINGS

TERACEL® Absorbable Haemostat is presented in its sterile original pack, cannot be re-sterilized and used. Unpacked, used or damaged products should not be used in any case. Use on contaminated wounds without cleaning can cause complications. Contaminated areas/wounds should definitely be cleaned prior to use of **TERACEL®** Absorbable Haemostat. Haemostatic properties **TERACEL®** Absorbable Haemostat gives the best results when applied dry. Therefore; it should not be damped nor moistened with water, salt solutions or other liquids prior to and during application. **TERACEL®** Absorbable Haemostat can be left on tissue but should be removed after bleeding stops. **TERACEL®** Absorbable Haemostat should be removed in any areas of use in or around bone foramen, bony areas, spinal cord or optical nerve or chiasm; no matter what type of surgical application is done. Haemostat can cause paralysis and/or damages on nerves as it swells and generates

pressure. **TERACEL®** Absorbable Haemostat should not be soaked in any kinds of materials such as antimicrobial agents or tampon/haemostatic substances. **TERACEL®** Absorbable Haemostat should not be used as an anti-adhesive substance. Sterility is preserved only when opened under sterile conditions. Do not re-sterilize.

Do not use opened or damaged products.

PRECAUTIONS

For homeostasis, **TERACEL®** Absorbable Haemostat should be tightly held on area of bleeding until bleeding stops and used in amount required. In order to ease absorption and minimize foreign substance reactions; the remaining amount should be removed before surgical area is closed. During urological operations; amount of material used should be at minimum levels and occlusion possibilities in urethra, ureter or catheter, which can occur as a result of dislocation of **TERACEL®** Absorbable Haemostat, should be watched out. Since absorption of **TERACEL®** Absorbable Haemostat can be chemically prevented in areas, where cautery is applied; it should not be used prior to application of silver nitrate or other scabbing chemicals. All measures should be taken during otorhinolaryngologic surgery; in order to prevent **TERACEL®** Absorbable Haemostat from being inhaled. **TERACEL®** Absorbable Haemostat should not be too tight when used for bandaging during vascular surgical operations. Discard used products in medical waste container.

ADVERSE REACTIONS

Liquid "Encapsulation", and foreign body reactions were reported. Besides; stentonic effect has been observed, where Oxidized Cellulose Haemostat is used as bandage during vascular surgery. Although a direct relation between use of **TERACEL®** Absorbable Haemostat and stenosis could not be found; high attention should still be paid and application as tight bandages should be avoided. Paralysis and nerve damage have been determined when Oxidized Cellulose Haemostat is used in and around bone foramina, bony boundary areas, spinal cord and/or optical nerve and chiasm. While many of aforementioned reports are related to laminectomy; paralysis occasions have been reported to be related also to other procedures. Case of blindness has been observed when Oxidized Cellulose Haemostat is placed in anterior cranial fossa for surgical repair of torn left frontal lobe. It has been reported that difficulty in passing of urine through urethra after prostatectomy and possibility of extension in drainage can occur in cholecystectomy operations (amputation of gallbladder). Also there is a report indicating that ureter is blocked after kidney resection (amputation of kidney), where post-operative catheterization is required.

Occasional reports of "burning" and "stinging" sensations and sneezing when **TERACEL®** Absorbable Haemostat has been used as packing in epistaxis, are believed to be due to the low pH of the product. Burning has been reported when **TERACEL®** Absorbable Haemostat was applied after nasal polyp removal and after haemorrhoidectomy. Headache, burning, stinging, and sneezing in epistaxis and other rhinological procedures, and stinging when **TERACEL®** Absorbable Haemostat was applied on surface wounds (varicose ulcerations, dermabrasions, and donor sites) also have been reported.

DOSSAGE AND ADMINISTRATION

TERACEL® Absorbable Haemostat is presented in sterilized package and sterile technique should be used while unpacking the product. Minimum amount should be put on bleeding area in appropriate measures and kept on tissue until bleeding stops. Amount of use depends on severity and nature of bleeding. It is known that bleeding preventing effect of **TERACEL®** Absorbable Haemostat is higher during dry application. Moistening the

material with water or salt solutions is not recommended.

STERILITY

TERACEL® Absorbable Haemostat sterilized by Gamma irradiation.

Do not re-sterilize!

Do not use if package is opened or damaged!

STORAGE

TERACEL® Absorbable Haemostat should be preserved in its original pack; in dark, dry places, and under controlled conditions (Between 5°C-25°C). Expire date of product is on its pack. **TERACEL®** Absorbable Haemostat should not be used after expire date.

HOW SUPPLIED

TERACEL® Absorbable Haemostat is presented as sterilized and in its original pack. It has different types of dimensions.

SYMBOLS USED ON THE LABELS

	CE Symbol and identification number of the notified body. The product satisfies the Essential Requirements of directive 93/42/EEC on medical devices.
	Manufacturer
	Date of manufacture (Year/Month)
	Use by expiration date (Year/Month)
	Lot Number
	Product Code
	R : Sterilization Method - Gamma Irradiation
	Caution
	Consult instructions for use
	Keep away from sunlight
	Protect from humidity
	Store between : 5°C - 25°C
	Single use
	Do not re-sterilize
	Do not use if package is damaged

KATSAN KATGUT SANAYI VE TİCARET A.Ş.
1571 Sk. No:22 35110 Cinarlı Izmir - TÜRKİYE
Tel: +90 (232) 486 0910 Fax: +90 (232) 486 3379
Web: www.katsanas.com E-mail: info@katsanas.com