## KATSAN ALCALACTONE PLUS®

Antibakteriyel Poliglolikolik asit (75%) -ko-kaprolaktone (25%) (PGCL)

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## KATSAN ALCALACTONE PLUS®

Antibacterial polyglycolic acid (75%)-co-caprolactone (25%) (PGCL)

Sterile absorbable antibacterial surgical suture

### DESCRIPTION

**KATSAN ALCALACTONE PLUS** sutures are synthetic absorbable sterile surgical sutures which are made of Polyglycolide and caprolactone. **KATSAN ALCALACTONE PLUS** sutures contain "Chlorhexidine diacetate" a broad-spectrum antibacterial agent.

**KATSAN ALCALACTONE PLUS** sutures when implanted into a living organism, it is absorbed by that organism and causes no undue tissue irritation. **KATSAN ALCALACTONE PLUS** sutures have been found to be non-antigenic, non-pyrogenic and elicit only a mild tissue reaction during absorption.

**KATSAN ALCALACTONE PLUS** sutures meet all requirements for synthetic absorbable surgical sutures specified at the European Pharmacopoeia (EP) and United States Pharmacopoeia (USP).

In order to ease of operation, **KATSAN ALCALACTONE PLUS** sutures are coloured D&G violet No: 2 or undyed (beige, natural). **KATSAN ALCALACTONE PLUS** sutures have been presented in USP 5/0 – USP 2 gauge sizes in a variety of lengths from 10 cm to 150 cm, as non-needled or attached to stainless steel needles of varying types and sizes.

### INDICATIONS

**KATSAN ALCALACTONE PLUS** sutures are indicated general ligation and/or soft tissue approximation which absorbable suture is needed.

### ACTIONS

**KATSAN ALCALACTONE PLUS** sutures elicit a minimal acute inflammatory reaction in tissue and ingrowth of fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of **KATSAN ALCALACTONE PLUS** suture occurs by means of hydrolysis, gradually and decreases the strength in the body. The absorption process begins a loss of tensile strength followed by a loss of mass. Implantation studies indicate that **KATSAN ALCALACTONE PLUS** suture retains approximately 60% of the original tensile strength at 1 week and in excess of 30%-40% at 2 weeks post implantation. Absorption of **KATSAN ALCALACTONE PLUS** suture is essentially complete between 90 to 120 days. Wound support period of **KATSAN ALCALACTONE PLUS** is approximately 21 days.

### CONTRAINDICATIONS

**KATSAN ALCALACTONE PLUS** sutures, being absorbable, are not be used where extended approximation of tissue is required.

### WARNINGS

Users should be familiar with surgical procedures and techniques covering absorbable sutures before employing **KATSAN ALCALACTONE PLUS** sutures for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material is used. Users should consider the in vivo performance, under "ACTION" section, when selecting a suture.

As with any foreign body, in surgery of the urinary or biliary tracts, attention should be taken to avoid prolonged contact of any suture with salt solutions, as may result in calculus formation.

As any foreign material in the presence of bacterial contamination may enhance bacterial infectivity, acceptable surgical practice should be followed with respect to drainage

and closure of contaminated or infected wounds.

The use of this suture may be inappropriate in elderly malnourished or debilitated patients, or in patients suffering from conditions which may delay wound healing.

**KATSAN ALCALACTONE PLUS** is an absorbable suture material, the use of supplemental nonabsorbable sutures should be considered by the surgeon in the closure of the sites which may undergo expansion, stretching, or distention, or which may require additional support.

The safety and effectiveness of **KATSAN ALCALACTONE PLUS** sutures have not been established in neural tissue, adult cardiovascular tissue for use in microsurgery. Under certain circumstances, notably orthopaedic procedures, immobilization by external support may be employed at the discretion of the surgeon.

Sterility is preserved only when opened under sterile conditions.

Do not re-sterilize.

Store at room temperature and avoid prolonged exposure to elevated temperatures.

Discard opened packages and unused sutures.

Do not use opened or damaged products.

### PRECAUTIONS

In handling **KATSAN ALCALACTONE PLUS** or any other suture material, care should be taken to avoid damage from handling. Especially avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

In order to adequate knot security requires the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstance and the experience of the surgeon.

Conjunctival and vaginal mucosal sutures remaining in place for extended periods may be associated with localized irritation. Subcuticular sutures should be placed as deeply as possible in order to minimize the erythema and induration normally associated with absorption.

Acceptable surgical practice should be followed with respect to drainage and closure of infected wounds.

In order to prevent damaging needle points and swage areas, the needle should be grasped on its flattened area which is one-third (1/3) to one-half (1/2) of the distance from the swaged area of the needle. Reshaping needle may result in loss of strength of the needles and may result in less resistance of the needle to bending and breaking. To prevent inadvertent needle sticks, users should pay attention when handling surgical needles. Discard used needles and PGCL Plus surgical suture in "sharps" container.

Do not use expired products.

### ADVERSE REACTIONS

Adverse effects associated with the use of this device include wound dehiscence, failure to provide adequate wound support in closure of the sites where expansion, stretching, or distention occur, failure to provide adequate wound support in elderly, malnourished or debilitated patients or in patients suffering from conditions which may delay wound healing, infection, minimal acute inflammatory tissue reaction, localized irritation when skin sutures are left in place for greater than 7 days, suture extrusion and delayed absorption in tissue with poor blood supply, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, and transitory local irritation at the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood borne pathogens.

### STERILITY

**KATSAN ALCALACTONE PLUS** sutures are sterilized by Ethylene Oxide.

Do not re-sterilize!

Do not use if package is opened or damaged!

Discard opened unused sutures.

### STORAGE

Recommended storage conditions: Store between 5°C-25°C, away from moisture and direct heat.

Do not use after expiry date!


Protect from humidity.

### HOW SUPPLIED

**KATSAN ALCALACTONE PLUS** sutures are available as sterile, monofilament dyed (violet) or undyed (beige, natural) strands in sizes USP 5/0 through USP 2 (metric sizes 1-5) , in a variety of lengths, with or without various needles.

Sutures are presented in one and two dozen boxes.

### SYMBOLS USED ON THE LABELS

	→	CE Symbol and identification number of the notified body. The product satisfies the Essential Requirements of directive 93/42/EEC on medical devices.
	→	Manufacturer
	→	Date of manufacture (Year/Month)
	→	Use by expiration date (Year/Month)
	→	Lot Number
	→	Product Code
	→	EO : Sterilization Method - Ethylene Oxide
	→	Caution
	→	Consult instructions for use
	→	Keep away from sunlight
	→	Protect from humidity
	→	Store between : 5°C - 25°C
	→	Single use
	→	Do not resterilize
	→	Do not use if package is damaged

 **KATSAN KATGUT SANAYI VE TICARET A.S.**  
1571 Sk. No:22 35110 Cinarlı Izmir - TÜRKİYE  
Tel: +90 (232) 486 0910 Fax: +90 (232) 486 3379  
Web: www.katsanas.com E-mail: info@katsanas.com