### TERAMESH®

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#### TANIM

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Partially Absorbable Mesh

#### DESCRIPTION

**TERAMESH®** is composed of absorbable Poly (glycolide-co-caprolactone, 75/25) monofilament fiber and polypropylene non-absorbable monofilament fibers. FDA-approved D&C violet No.2 is used in coloured fibers.

#### INDICATIONS

**TERAMESH®** is indicated in hernia, eventration, rectal and genital urinary prolapses cases; also used for treatment of abdominal and incisional hernias and suitable for use in open or laparoscopic techniques.

#### ACTIONS

**TERAMESH®** is used for abdominal wall stabilization in cases such as hernia repair. The absorbable fiber part of material ensures that the propylene part remains stabilized and provides ease of usage. **TERAMESH®**, which is considered as partially absorbable when compared to completely polypropylene meshes, maintains strongly effective advantages on patient's comfort.

#### CONTRAINDICATIONS

**TERAMESH®** should be separated from peritoneum and abdominal cavity and not be used for people allergic to polypropylene. **TERAMESH®** is not appropriate for being placed in groin canal in straight mesh form as stopper. **TERAMESH®** should not be used after gastrointestinal system is opened accidentally or on purpose during operation. In such cases; it may cause contamination of mesh and eventually an infection that can require mesh to be removed.

#### WARNINGS

**TERAMESH®** is presented in its sterile original pack and should not be re-sterilized. Unpacked or damaged products should not be used in any case. In cases that **TERAMESH®** is used for babies, kids, pregnant or pregnancy-planning women; surgeons should be aware that the product

does not have elasticity enough to expand as much as patient's body growths. Acute or chronic; local or systemic infections; muscle based, neurological or vascular incompetence; any kind of natural diseases can effect function of implant. Some patients can face foreign substance reactions because of its use. **TERAMESH®** should not be used in contaminated surgical areas. In such conditions; strong infections, which can require material to be removed out, can be observed. Prior to use of product; surgeon should have an adequate knowledge about techniques of usage and application for non-absorbable or partially absorbable meshes. Sterile technique should be used while unpacking the product and not be unpacked until required for use. The product cannot be re-used nor re-sterilized for second use. Sterility is preserved only when opened under sterile conditions.

Do not re-sterilize.  
Do not use opened or damaged products.

#### ADVERSE REACTIONS

Possible adverse reactions related to **TERAMESH®** implantation are the ones seen generally on surgically implanted materials and consist of inflammation, seroma formation, adhesion formation, fistula formation, extrusion and infection to become intensified.

#### STERILITY

**TERAMESH®** is sterilized by Ethylene Oxide. Do not re-sterilize!  
Do not use if package is opened or damaged!

#### STORAGE

**TERAMESH®** should be preserved in its original pack; in dark, dry places, and under controlled conditions (Between 5°C-25°C). Expire date of product is on its pack. **TERAMESH®** should not be used after expire date. Discard used products in medical waste container.

#### HOW SUPPLIED

**TERAMESH®** is presented as sterilized and in its original pack. It has different types of dimensions.

#### SYMBOLS USED ON THE LABELS

	→	CE Symbol and identification number of the notified body. The product satisfies the Essential Requirements of directive 93/42/EEC on medical devices.
	→	Manufacturer
	→	Date of manufacture (Year/Month)
	→	Use by expiration date (Year/Month)
	→	Lot Number
	→	Product Code
	→	EO : Sterilization Method - Ethylene Oxide
	→	Caution
	→	Consult instructions for use
	→	Keep away from sunlight
	→	Protect from humidity
	→	Store between : 5°C - 25°C
	→	Single use
	→	Do not re-sterilize
	→	Do not use if package is damaged

**KATSAN KATGUT SANAYI VE TICARET A.S.**  
1571 Sk. No:22 35110 Cinarli Izmir - TURKIYE  
Tel: +90 (232) 486 0910 Fax: +90 (232) 486 3379  
Web: www.katsanas.com E-mail: info@katsanas.com